


Gas Specific Quick Connect



FEATURES

- Use to mount Medical Gas Pressure/Vacuum Switches and Gauges
- Allows Alarm Testing without Service Disruption
- Simplifies Servicing Pressure Switches and Gauges

SPECIFICATIONS

The gas specific demand check fitting (as required by NFPA 99) shall be designed to allow quick removal of pressure/vacuum switches and gauges from the pipeline. Facilitates replacement and ease of periodic testing without disrupting the system.

The gas specific demand check assembly shall consist of 1/4" chrome plated brass female check unit and a mating chrome plated brass fitting with 1/8" NPT thread.

PART NUMBERS

<input type="checkbox"/> 261977	Evacuation	<input type="checkbox"/> 261924	N2**
<input type="checkbox"/> 261976	CO2	<input type="checkbox"/> 261932	Vacuum
<input type="checkbox"/> 261922	Air	<input type="checkbox"/> 261933	Oxygen
<input type="checkbox"/> 261923	N2O		

**Nitrogen (N2) requires two pressure switches. Replace the tee with a cross to install two pressure switches

TRUSTED BRANDS OF OHIO MEDICAL



Ohio Medical Corporation - 1111 Lakeside Drive - Gurnee, IL 60031-4099 USA
 TollFree: 800-448-0770 - Phone: 847-855-0500 - Fax: 847-855-6300
www.ohiomedical.com